




Wieslab® Product Range

Leading Specialised Autoimmune Tests

The Wieslab® line of diagnostic products has been developed with leading clinicians in the field. This ensures the kits always give optimal performance and is the reason why the Wieslab® kits are used by the top laboratories around the world.

- » ANCA assays calibrated against international standards
- » Assessment of complement function
- » Unique products for diagnosis of rare diseases
- » Vasculitis diagnosis – capture format for improved sensitivity/specificity



Analyte Antibody class detected Product code	Plate coating	Calibrators	Units	Quant	Qual	Incubation time	Tests	CE
ANCA panel kit IgG PAN106	AZ, BPI, CH, EL, LF, LY	Ratio	N/A	No	Yes	30+30+60 min	96	x
ANCA screen IgG CP111	MPO, PR3	Ratio	N/A	No	Yes	10+10+20 min	96	x
Vasculitis screen IgG GCP-CAP	GBM, MPO, PR3	Ratio	N/A	No	Yes	10+10+20 min	96	x
Anti-GBM IgG GP104X	GBM (bovine α3 chain)	Yes	10-320 (U/mL)	Yes	No	30+30+60 min	96	x
Anti-GBM, ANCA screen IgG GCP100	GBM, MPO, PR3	Ratio	N/A	No	Yes	10+10+20 min	96	x
Capture MPO ANCA IU IgG Cap MPO IU	MAB + human MPO	Yes	0-200 (IU/mL) ¹	Yes	No	60+30+30 min	96	x
Capture PR3 ANCA IU IgG Cap PR3 IU	MAB + human PR3	Yes	0-200 (IU/mL) ¹	Yes	No	60+30+30 min	96	x
MPO ANCA IU IgG MPO IU	Human MPO	Yes	0-200 (IU/mL) ¹	Yes	No	60+30+30 min	96	x
PR3 ANCA IU IgG PR3 IU	Human PR3	Yes	0-200 (IU/mL) ¹	Yes	No	60+30+30 min	96	x
ASCA IgA ASCA150	Mannan	Yes	10-160 (U/mL)	Yes	Yes	30+30+60 min	96	x
ASCA IgM ASCA151	Mannan	Yes	10-160 (U/mL)	Yes	Yes	30+30+60 min	96	x
Complement system screen COMPL300	Human IgM + LPS + mannan	2 points	0-100%	No	Yes	60+30+30 min	96	x
Complement system AP COMPLAP330	LPS	2 points	0-100%	No	Yes	60+30+30 min	96	x
Complement system CP COMPLCP310	Human IgM	2 points	0-100%	No	Yes	60+30+30 min	96	x
Complement system MP COMPLMP320	Mannan	2 points	0-100%	No	Yes	60+30+30 min	96	x
SS-A p200 IgG SSAp200	Recombinant peptide	Ratio	N/A	No	Yes	60+30+30 min	96	x
Alport's syndrome ² N/A ALP105	N/A	N/A	N/A	N/A	Yes	See IFU	N/A	x
hCOMP IgG COMP200	Human cartilage (COMP)	Yes	10-80 ng/mL	Yes	No	See IFU	96	
sGAG ³ N/A GAG201	N/A	Yes	12.5-400 µg/mL	Yes	No	See IFU	200	

¹ Calibrated against International Standards (IU).

² Immunohistochemical staining kit.

³ Dye-binding assay.

Wieslab® product range • Leading specialised autoimmune tests



Novo Healthcare Ltd.

Corporate Office: Global House, Tower 1, 100 Broad Street, 4th Floor, London, E14 4EU, UK
 11th Floor, Central House, Broad Street, London, E14 4EU, UK
 Tel: +44 (0)20 7464 7222 | Fax: +44 (0)20 7464 7271 | Email: info@novohc.com

EuroDiagnostica